


STI Specimen Collection Requirements

January 2024



- Patient label MUST not cover blue label on tube.
- Leave part of fill window visible.
- Urine samples MUST be filled between MAX-MIN lines on fill window. Do not overfill or underfill, samples will be rejected.
- Ensure approved source is entered electronically when order is placed.

Specimen Type	CTPCR	NGPCR	TVPCR	STI2	STI3
CERVICAL	✓	✓	✓	✓	✓
VAGINAL	✓	✓	✓	✓	✓
URINE	✓	✓	✓	✓	✓
THINPREP	✓	✓	✓	✓	✓
RECTAL	✓	✓		✓	
THROAT	✓	✓		✓	