







For any specimen collected using a container/swab outside of this list please call the NoCo Microbiology laboratory at 970-820-4156 for approval

Aerobic Swab - eSwab	
<p>Rapid Group A Strep (throat) by PCR</p> <p>Aerobic Culture Anaerobic Culture Burn Culture Ear Culture Eye Culture Genital Culture GC Culture MRSA Screen Strep Screen Culture VRE Screen</p>	
BD NP swab w/3 ml UTM	
<p>Rapid Influenza A&B/RSV by PCR</p> <p>Expanded Respiratory Virus Panel by PCR</p>	
COPAN Swab	
<p>Prenatal Group B Strep by Amplification</p>	
Minitip Swab	
<p>Nasopharyngeal Culture Bordetella Pertussis by Amplification</p>	

Once printed this document may no longer be current.
Please refer to the online Laboratory Services Directory (JDOS) for current information.

<p>BD Affirm Collection Device</p> <p>Vaginitis testing (automated wet prep) screens for: Candida, Trichomonas, Gardnerella</p>	
<p>Gen-Probe Aptima</p> <p>Urine Specimen Male and Female GC / Chlamydia</p>	
<p>Gen-Probe Aptima Unisex</p> <p>Endocervical Male Urethral GC / Chlamydia (Can also be used for alternative site testing for GC/Chlamydia)</p>	
<p>Para - Pak</p> <p>Ova and Parasite Exam Giardia and Cryptosporidium Antigen</p> <p>*DO NOT Refrigerate</p>	
<p>Para – Pak C&S</p> <p>Stool Pathogens by PCR Giardia and Cryptosporidium Antigen</p> <p>*REFRIGERATE</p>	

Once printed this document may no longer be current.
Please refer to the online Laboratory Services Directory (JDOS) for current information.

Sterile Container		
<p>Clostridium Difficile Urine Fecal WBCs</p> <p>*REFRIGERATE</p>		