

Light Blue Top (Plastic-Citrated) Tube Draw Volume Guideline

BD Vacutainer® Plus Plastic Citrate Tube

BD Vacutainer® Plus Plastic Citrate Tube Draw Volume Guide

*Sufficient volume achieved
if blood drawn falls above
minimum fill indicator.
For blood transfer, **do not**
fill above illustrated dashed
maximum line.*

Note: The quantity of blood drawn
into evacuated tubes varies with
altitude, ambient temperature,
barometric pressure, tube age,
venous pressure and filling technique.

*According to CLSI guideline,
Dec. 2003, Doc. H1-A5, Vol. 23, No. 33.

