

KIDNEY TRANSPLANT Requisition for HLA Testing

Patient Name: _____ Ordering Physician: _____

Date of Birth: _____ Race: _____ Sex: _____ ABO: _____

WVUM MRN: _____ SSN: _____

Diagnosis/Code: _____ Transplant type: **KIDNEY**

Collection Date: _____ Collection Time: _____ Phlebotomist: _____

— INITIAL PATIENT EVALUATION (WVUH Epic # LAB123997)

1 RED TOP CLOT TUBE; 10 cc (HLA Antibody testing)

4 ACD YELLOW TOPS; 40 cc (HLA typing)

Note: This includes: HLA Antibody testing Luminex PRA (Class I and Class II) (LAB123994), HLA Antibody testing Luminex SAB (Class I and Class II) (High definition) (LAB123995), and HLA Typing Class I and Class II (LAB123996)

— LUMINEX MONTHLY SERUM SCREENING (PRA) / SERUM STORAGE (WVUH Epic # LAB123994)

1 RED TOP CLOT TUBE; 10 cc (HLA Antibody testing / Crossmatching)

— HLA TYPING CLASS I & II (WVUH Epic # LAB123996)

4 ACD YELLOW TOPS; 40 cc (HLA typing)

— LIVING DONOR WORKUP—RECIPIENT (INITIAL) (WVUH Epic # LAB123991)

1 RED TOP CLOT TUBE; 10 cc (HLA Antibody testing / Crossmatching)

4 ACD YELLOW TOPS; 40 cc (if recipient **NOT** previously HLA typed)

— LIVING DONOR WORKUP—POTENTIAL DONOR (INITIAL OR FINAL) (WVUH Epic # LAB123797)

4 ACD YELLOW TOPS; 40 cc *Note: testing to be ordered on donor and not recipient*

Recipient Name _____ DOB: _____

Relationship: _____

— LIVING DONOR WORKUP—RECIPIENT (FINAL) (WVUH Epic # LAB123992)

1 RED TOP CLOT TUBE; 10 cc (HLA Antibody testing / Crossmatching)

— LUMINEX HIGH DEFINITION (SAB): CLASS I __ CLASS II __ (WVUH Epic # LAB123995)

1 RED TOP CLOT TUBE; 10 cc (HLA Antibody testing / Crossmatching)

— LUMINEX DSA TESTING—post-transplant patients only (WVUH Epic # LAB123993)

1 RED TOP CLOT TUBE; 10 cc (HLA Antibody testing)

— DECEASED (CADAVER) DONOR CROSSMATCH (WVUH Epic # LAB123798)

1 RED TOP CLOT TUBE; 10 cc (HLA Antibody testing)

LAB STAFF:

Do not centrifuge

Store at Room Temperature

Ship via Overnight Fed Ex

Please send completed requisition and blood samples to:

a. US Mail: Use UN3373 Diamond Label or mark outer package "DIAGNOSTIC SPECIMEN"

b. Overnight: UPS or Federal Express – label outer package with UN3373 Diamond Label