

# NORTH MISSISSIPPI MEDICAL CENTER

Tupelo, Mississippi

## Bone Marrow Requisition

PATIENT IDENTIFICATION

### PATIENT INFORMATION:

Name: \_\_\_\_\_ Doctor: \_\_\_\_\_  
( ) Inpatient Room #: \_\_\_\_\_ MR#: \_\_\_\_\_ Date: \_\_\_\_\_  
( ) Hematology/Oncology Clinic

SITE: [ ] Right Iliac Crest [ ] Left Iliac Crest [ ] Other:

### WORKING CLINICAL DIAGNOSIS (Reason for bone marrow exam):

- [ ] Anemia [ ] Rule out Granuloma [ ] Thrombocytopenia  
[ ] Pancytopenia [ ] Possible Infection [ ] Acute Leukemia  
[ ] Possible Lymphoma [ ] Lymphoma Staging  
[ ] Possible Metastatic Tumor (Tumor Type: \_\_\_\_\_)  
[ ] Patient Immunosuppressed  
[ ] Rule out Myeloproliferative Disorder  
[ ] Post-Induction Chemotherapy for Acute Leukemia  
[ ] Monoclonal Protein ( ) Serum ( ) Urine  
[ ] Previous Chemotherapy (Give date of therapy: \_\_\_\_\_)  
[ ] Other (Detailed)

### MATERIAL:

- [ ] Marrow Smears (Number: \_\_\_\_\_)  
[ ] Aspirate Material for Clot Preparation  
[ ] Touch Preps from Biopsy (Number: \_\_\_\_\_)  
[ ] Bone Marrow Biopsy  
[ ] Bone Marrow (0.5-1.0 mls) in an ACD (yellow stoppered) Tube  
[ ] Send for Flow Cytometry (Quest-35080)(SQ-QLCMS)  
[ ] Bone Marrow (1 ml) in Sodium Heparin (green stoppered) Tube  
[ ] Send for:  
[ ] Cytogenetics - (Quest-14600X)(SQ-QCHHB)  
[ ] FISH Panel - Myelodysplasia (Quest-19799X)(SQ-QFMDS)  
[ ] FISH Panel - Acute Myeloid Leukemia (Quest-90871)(SQ-QFAML)  
[ ] FISH Panel - AML M3,PML/RARA, Translocation 15/17 (Quest 14617) (SQ-Order Flexitest)  
[ ] FISH Panel - Acute Lymphoid Leukemia (Quest-40052)(SQ-QFBAL)  
[ ] FISH Panel - Plasma Cell Myeloma (Quest-92497)(SQ-QFPCP)  
[ ] FISH Panel - High Grade B-Cell Lymphoma (Quest-91606)(SQ-QFBLP)  
[ ] FISH Panel - Chronic Myelogenous Leukemia (9;22 translocation)(Quest-12070X)(SQ-QFISH)  
[ ] FISH Panel - Chronic Lymphocytic Leukemia (Quest-16864) (SQ-FCLL)  
[ ] Other Special Tests Desired:

Ordering Physician: \_\_\_\_\_ Date: \_\_\_\_\_ Time: \_\_\_\_\_

### (Attending Pathologist to modify tests above as deemed necessary)

[ ] 500 CELL DIFF [ ] 200 CELL DIFF [ ] NO DIFF



Lab Use Only	Stained	Unstained
Bone Marrow Smears		
Marrow Smear Reviewed by: _____, Pathologist		